

DRIVING INNOVATION. DELIVERING THROUGHPUT.

From the people who set the standard for Low Temperature Sterilization.

The ***New Standard***
for Delivering to the O.R.

V-PRO[®]maX²
Low Temperature Sterilization System



MAXIMUM THROUGHPUT

The V-PRO maX 2 Sterilizer delivers greater productivity for day-to-day OR needs than any sterilizer on the market.

Max Speed	For Non Lumen Devices	Max Capacity	For Non Lumen Devices
V-PRO^{maX} 2 Low Temperature Sterilization System V-PRO maX 2 Sterilizer Fast Non Lumen Cycle	 16 minute cycle time	V-PRO^{maX} 2 Low Temperature Sterilization System V-PRO maX 2 Sterilizer Non Lumen Cycle	 28 minute cycle time
Brand X STERRAD 100NX [®] Sterilizer Express Cycle	 24 minute cycle time	Brand X STERRAD 100NX [®] Sterilizer Standard Cycle	 47 minute cycle time

Popular devices validated

Experience the speed of the innovative Fast Non Lumen Cycle and process devices in just 16 minutes, such as


- Cameras
- Telescopes
- Batteries
- Laryngoscope blades



Full list of validated devices: steris.com/devicematrix


Max Ease

Simple, intuitive cycle selection via large icons make training and day-to-day use easy.



Minimize Aborted Cycles

Delayed Start and Moisture Check features save time and money.



Hands-Free Loading

Innovative foot pedal makes loading quick, safe and easy.

Less Sterilant, More Value

Fast Non Lumen Cycles uses half the sterilant compared to all other V-PRO Sterilizer cycles, lowering your cost-per-cycle.



Know You're Running Low

RFID technology notifies you of remaining sterilant and cup expiration to help reduce waste.

PRO-LITE[™] Sterilization Trays

Designed specifically for V-PRO Sterilizers, optimize loading capacity and increase weight capacity.



INCREASED PRODUCTIVITY.
REDUCED EXPENSE.

* Savings based on U.S. Government contract pricing for STERRAD 100NX in 2014 over a 7 year lifetime.
** With Celerity 20 HP Biological Indicator



35% less
consumable
cost*



27%
less service
cost*



20%
faster
BI read**

V-PRO maX 2 Sterilizer + Celerity™ 20 HP Biological Indicator = MAXIMUM THROUGHPUT — **Fastest Results** in the Industry



celerity™ 20 HP Biological Indicator

- 20% faster than leading competitor
- Accelerated load release means improved productivity
- Innovative glass-free design increases staff safety

COMPLIANCE MADE SIMPLE



Driven by **real-time data** delivered from your V-PRO maX 2 Sterilizer and Celerity HP Incubator to your facility's tracking system, **ConnectAssure** Technology allows your staff to focus on what matters most — **supporting the OR.**



- Complete sterilization records in one place
- Reminders for regular maintenance tasks
- Essential reporting for compliance
 - Demonstrate compliance over a defined period
 - Eliminate lost or incomplete sterilization records

No tracking system? No problem. ConnectAssure Technology has a robust interface that allows all the benefits of the software without a facility tracking system.

To learn more about V-PRO Sterilizers, contact your
STERIS Representative or visit STERIS.com/VPROmaX2



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com